


Product Datasheet

Tau MAPT Rabbit Monoclonal Antibody (orb547051)

Description	Tau MAPT Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, WB
Immunogen	A synthesized peptide derived from human Tau
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	78928 MW
Uniprot ID	P10636
Dilution Range	WB 1:1000-1:2000ICC/IF 1:50-1:200FC 1:100
Expiration Date	12 months from date of receipt.



WB analysis
of Tau
expression
in Mouse
r...