



Product Datasheet

KRAS+HRAS+NRAS Rabbit Monoclonal Antibody (orb547055)


Description	KRAS+HRAS+NRAS Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IP, WB
Immunogen	A synthesized peptide derived from human KRAS+HRAS+NRAS
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	21229 MW
Uniprot ID	P01116, P01111, P01112
Dilution Range	WB 1:500-1:2000 ICC/IF 1:50-1:200 IP 1:50 FC 1:200
Expiration Date	12 months from date of receipt.



KRAS+HRAS+NRAS was immunoprecipitated fr...



WB analysis of KRAS+HRAS+NRAS expression...



Immunofluorescent analysis of HeLa cells...