


Product Datasheet

Ras NRAS Rabbit Monoclonal Antibody (orb547056)

Description	Ras NRAS Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	A synthesized peptide derived from human Ras
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	21229 MW
Uniprot ID	P01116, P01111, P01112
Dilution Range	WB 1:10000-1:20000IP 1:20FC 1:50
Expiration Date	12 months from date of receipt.



WB analysis
of Ras
expression
in (1) C6 ...