


Product Datasheet

Huntingtin Monoclonal Antibody (orb547087)

Description	Huntingtin Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human Huntingtin May play a role in microtubule-mediated transport or vesicle function.
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P42858
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:100-1:500FC 1:100
Expiration Date	12 months from date of receipt.



WB analysis
of HTT
expression
in HeLa Ce...