


## Product Datasheet

# alpha + beta Synuclein SNCA Rabbit Monoclonal Antibody (orb547159)

|                            |                                                                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | alpha + beta Synuclein SNCA Rabbit Monoclonal Antibody                                                                                     |
| <b>Species/Host</b>        | Rabbit                                                                                                                                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                                                                               |
| <b>Tested Applications</b> | FC, ICC, IF, WB                                                                                                                            |
| <b>Immunogen</b>           | A synthesized peptide derived from human alpha + beta Synuclein                                                                            |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                     |
| <b>Concentration</b>       | Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.                                               |
| <b>Storage</b>             | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                      |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                 |
| <b>Clonality</b>           | Monoclonal                                                                                                                                 |
| <b>MW</b>                  | 14460 MW                                                                                                                                   |
| <b>Uniprot ID</b>          | <b>Q16143</b> , <b>P37840</b>                                                                                                              |
| <b>Dilution Range</b>      | WB 1:1000-1:2000ICC/IF 1:50-1:200FC 1:50                                                                                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                            |



WB  
analysis of  
alpha +  
beta  
Synuclein  
ex...