


Product Datasheet

TNFAIP3/A20 Rabbit Monoclonal Antibody (orb547168)

Description	TNFAIP3/A20 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human TNFAIP3
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	89614 MW
Uniprot ID	P21580
Dilution Range	WB 1:1000-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200FC 1:50
Expiration Date	12 months from date of receipt.



WB
analysis of
TNFAIP3
expression
in Jur...