



Product Datasheet

Hemoglobin subunit alpha HBA1 Rabbit Monoclonal Antibody (orb547176)


Description	Hemoglobin subunit alpha HBA1 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP, WB
Immunogen	A synthesized peptide derived from human Hemoglobin subunit alpha
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	15258 MW
Uniprot ID	P69905
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200IP 1:50FC 1:50
Expiration Date	12 months from date of receipt.



WB analysis of Hemoglobin subunit alpha ...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...