


Product Datasheet

TBP/Tata Binding Protein Tbp Rabbit Monoclonal Antibody (orb547226)

Description	TBP/Tata Binding Protein Tbp Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human TBP
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	37698 MW
Uniprot ID	P20226
Dilution Range	WB 1:1000-1:2000
Expiration Date	12 months from date of receipt.



WB analysis of TBP expression in K562 Ce...