



Product Datasheet

iNOS NOS2 Rabbit Monoclonal Antibody (orb547276)

Description	iNOS NOS2 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	A synthesized peptide derived from human iNOS
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	131117 MW
Uniprot ID	P35228
Dilution Range	WB 1:500-1:1000 ICC/IF 1:50-1:200 IP 1:50
Expiration Date	12 months from date of receipt.



WB analysis of
iNOS expression in
Human ...



Immunofluorescent
analysis of
Raw264.7 c...