



Product Datasheet

ATF4 Rabbit Monoclonal Antibody (orb547278)


Description	ATF4 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human ATF4
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	38590 MW
Uniprot ID	P18848
Dilution Range	WB 1:500-1:1000IHC 1:50-1:200ICC/IF 1:100-1:500IP 1:50FC 1:50
Expiration Date	12 months from date of receipt.



WB analysis of ATF4 expression in (1) He...



Immunofluorescent analysis of 293 cells,...



Immunohistochemical analysis of paraffin...