


Product Datasheet

AGTR2/Angiotensin II Type 2 Receptor Rabbit Monoclonal Antibody (orb547313)

Description	AGTR2/Angiotensin II Type 2 Receptor Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthesized peptide derived from human AGTR2
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	41184 MW
Uniprot ID	P50052
Dilution Range	WB 1:1000-1:2000IP 1:50
Expiration Date	12 months from date of receipt.



WB
analysis of
AGTR2
expression
in HepG2...