


Product Datasheet

ErbB (HER3) Rabbit Monoclonal Antibody (orb547362)

Description	ErbB (HER3) Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthesized peptide derived from human ErbB (HER3)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	148098 MW
Uniprot ID	P21860
Sensitivity	10pg/ml
Dilution Range	WB 1:500-1:2000IP 1:30
Expiration Date	12 months from date of receipt.



WB analysis
of ErbB
(HER3)
expression
in...