


Product Datasheet

CDK7 Rabbit Monoclonal Antibody (orb547374)

Description	CDK7 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human CDK7
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	39038 MW
Uniprot ID	P50613
Dilution Range	WB 1:1000-1:2000
Expiration Date	12 months from date of receipt.



WB analysis of CDK7 expression in MCF-7 ...