



Product Datasheet

TrkA+B+C NTRK1 Rabbit Monoclonal Antibody (orb547441)


Description	TrkA+B+C NTRK1 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human TrkA+B+C
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	87497 MW
Uniprot ID	Q16288 , P04629 , Q16620
Dilution Range	WB 1:1000-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200
Expiration Date	12 months from date of receipt.



WB analysis of TrkA+B+C expression in Hu...



Immunofluorescent analysis of Neuro-2a C...



Immunohistochemical analysis of paraffin...