



Product Datasheet

APOBEC3G Antibody (monoclonal, 6C2) (orb547443)


Description	APOBEC3G Antibody (monoclonal, 6C2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human APOBEC3G recombinant protein (Position: E191-N384).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	46 kDa
Uniprot ID	Q9HC16
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Immunocytochemistry/Immunofluorescence, 2µg/ml Flow Cytometry, 1-3µg/1x10 ⁶ cells
Expiration Date	12 months from date of receipt.



Flow
Cytometry
analysis of
THP-1 cells
u...



WB analysis
of APOBEC3G
using anti-
APOBE...



IF analysis of
APOBEC3G
using anti-
APOBE...