



Product Datasheet

Collagen III/COL3A1 Antibody (monoclonal, 9H9) (orb547486)

Description	Collagen III/COL3A1 Antibody (monoclonal, 9H9)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	E. coli-derived human Collagen III/COL3A1 recombinant protein (Position: D1222- E1455).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	180-190 kDa
Uniprot ID	P02461
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml
Expiration Date	12 months from date of receipt.



WB analysis of Collagen III/COL3A1 using...



IHC analysis of Collagen III/COL3A1 usin...



IHC analysis of Collagen III/COL3A1 usin...