



Product Datasheet

Nestin Rabbit Monoclonal Antibody (orb547493)

Description	Nestin Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	A synthesized peptide derived from human Nestin
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	177439 MW
Uniprot ID	P48681
Dilution Range	IHC 1:50-1:200 ICC/IF 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunofluorescent analysis of HeLa cells...