



Product Datasheet

S100B/S100 Beta Rabbit Monoclonal Antibody (orb547555)


Description	S100B/S100 Beta Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human S100B
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	10713 MW
Uniprot ID	P04271
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunofluorescent analysis of A375 cells...



WB analysis of S100B expression in A375 ...