


Product Datasheet

Mouse Efemp1 protein (orb54758)

Description	Full Length of Mature Protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal Flag-Myc-tagged
Note	For research use only
Application notes	N-terminal 10xHis-tagged and C-terminal Myc-tagged: N-terminal 6xHis-tagged1-518AA: 18-493AAFull Length : Full Length
Protein Sequence	QYTEETITYTQCTDGYEWDPIRQQCKDIDECDIVPDACKGGMKCVNHYGGYLCL PKTAQIIVNNEHPQQETPAAEASSGATTGTVAARSMATSGVVPGGGFMASATAV AGPEVQTGRNNFVIRRNAPDQRIIPSNPSHRIQCAAGYEQSEHNVCQDIDECTS GTHNCRTDQVCINLRGSFTQCCLPGYQKRGEQCVDIDECTVPPYCHQRCVNT GSFYCQCSPGFQLAANNYTCVDINECDASNQCAQQCYNILGSFICQCNQGYELS SDRLNCEDIDECRTSSYLCQYQCVNEPGKFSCMCPQGYEVVRSRTCQDINECET TNECREDEM CWNYHGGFRCPYRNPQCQDHYVLTSENRCVCPVSNMTCRELPQSI VYKYMSIRSDRSVPSDIFIQIATMIYANTINTFRIKSGNENGIFYLRQTSPVSAMLV LVKSLSGPREYIVDLEMLTVSSIGTFR TSSVLR LTIIVGPFSF
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Mammalian cell
MW	57.5 kDa
Uniprot ID	Q8BPB5
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Mouse
Efemp1
protei...