


Product Datasheet

Endothelin B Receptor EDNRB Rabbit Monoclonal Antibody (orb547584)

Description	Endothelin B Receptor EDNRB Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthesized peptide derived from human Endothelin B Receptor
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	49644 MW
Uniprot ID	P24530
Dilution Range	WB 1:500-1:2000IP 1:50
Expiration Date	12 months from date of receipt.



WB
analysis of
Endothelin
B Receptor
exp...