


Product Datasheet

PI3 Kinase p110 beta PIK3CB Rabbit Monoclonal Antibody (orb547605)

Description	PI3 Kinase p110 beta PIK3CB Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	A synthesized peptide derived from human PI3 Kinase p110 beta
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	122762 MW
Uniprot ID	P42338
Dilution Range	WB 1:500-1:2000IP 1:50FC 1:50
Expiration Date	12 months from date of receipt.



WB
analysis of
PI3 Kinase
p110 beta
expr...