



Product Datasheet

Hsp90 alpha + beta HSP90AA1 Rabbit Monoclonal Antibody (orb547614)


Description	Hsp90 alpha + beta HSP90AA1 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human Hsp90 alpha + beta
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	84660 MW
Uniprot ID	P08238, P07900
Dilution Range	WB 1:5000-1:10000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:50FC 1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...