


Product Datasheet

Troponin T TNNT2 Rabbit Monoclonal Antibody (orb547634)

Description	Troponin T TNNT2 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthesized peptide derived from human Troponin T
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	35924 MW
Uniprot ID	P45379
Dilution Range	WB 1:5000-1:20000IHC 1:50-1:200
Expiration Date	12 months from date of receipt.



WB
analysis of
Troponin T
expression
in ...