



Product Datasheet

NGFR/P75 Rabbit Monoclonal Antibody (orb547646)

Description	NGFR/P75 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human p75 NGF Receptor
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	45183 MW
Uniprot ID	P08138
Sensitivity	10pg/ml
Dilution Range	WB 1:10000-1:20000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



WB analysis of NGFR expression in (1) C6...