


Product Datasheet

Heparanase 1 HPSE Monoclonal Antibody (orb547701)

Description	Heparanase 1 HPSE Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthesized peptide derived from human Heparanase 1.
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q9Y251
Dilution Range	WB 1:500-1:2000IP 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of Heparanase 1 ex...