


Product Datasheet

Torsin A TOR1A Rabbit Monoclonal Antibody (orb547726)

Description	Torsin A TOR1A Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthesized peptide derived from human Torsin A
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	37809 MW
Uniprot ID	O14656
Sensitivity	> 5000 cells
Dilution Range	WB 1:500-1:2000FC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of Torsin A expres...