



Product Datasheet

Iba1 AIF1 Monoclonal Antibody (orb547733)


Description	Iba1 AIF1 Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human Iba1
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P55008
Dilution Range	WB 1:500-1:1000IHC 1:100-1:500ICC/IF 1:50-1:200IP 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of Iba1 expression...



Immunohistochemical analysis of paraffin...



Immunofluorescence staining of paraffin...