



Product Datasheet

beta III Tubulin TUBB3 Rabbit Monoclonal Antibody (orb547882)

Description	beta III Tubulin TUBB3 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human beta III Tubulin
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	50433 MW
Uniprot ID	Q13509
Dilution Range	WB 1:5000-1:10000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:50FC 1:50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of beta III Tubuli...