



Product Datasheet

GAP43 Rabbit Monoclonal Antibody (orb547884)

Description	GAP43 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human GAP43
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	24803 MW
Uniprot ID	P17677
Dilution Range	WB 1:10000-1:50000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:30FC 1:20
Expiration Date	12 months from date of receipt.



Western blot analysis of GAP43 expressio...



Immunohistochemical analysis of paraffin...