


Product Datasheet

TRP1 Monoclonal Antibody (orb547914)

Description	TRP1 Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	A synthesized peptide derived from human TRP1 Oxidation of 5,6-dihydroxyindole-2-carboxylic acid (DHICA) into indole-5,6-quinone-2-carboxylic acid. May regulate or influence the type of melanin synthesized.
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P17643
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200IP 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of TYRP1 expression...