


Product Datasheet

JAM1 Rabbit Monoclonal Antibody (orb547943)

Description	JAM1 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthesized peptide derived from human JAM1
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	32583 MW
Uniprot ID	Q9Y624
Sensitivity	10pg/ml
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of JAM1 expression...