

Product Datasheet

PRP19 PRPF19 Rabbit Monoclonal Antibody (orb548042)

Description PRP19 PRPF19 Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, ICC, IF, IHC, WB

Immunogen A synthesized peptide derived from human PRP19

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Note For research use only

Isotype Rabbit IgG


Clonality Monoclonal

MW 55181 MW


Uniprot ID **Q9UMS4**

Dilution Range WB 1:500-1:2000 IHC 1:50-1:200 ICC/IF 1:50-1:200 FC 1:50

Expiration Date 12 months from date of receipt.



Western blot analysis of PRP19 expression...



Immunohistochemical analysis of paraffin...