



Product Datasheet

PPT1 Antibody (monoclonal, 10F3) (orb548087)


Description	PPT1 Antibody (monoclonal, 10F3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human PPT1 (191-224aa KEDVYRNHSIFLADINQERGINESYKKNLMALKK), different from the related mouse and rat sequences by four amino acids.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	34 kDa
Uniprot ID	P50897
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Flow Cytometry, 1-3µg/1x10 ⁶ cells
Expiration Date	12 months from date of receipt.



Western blot analysis of PPT1 using anti...



IHC analysis of PPT1 using anti PPT1 ant...



IHC analysis of PPT1 using anti PPT1 ant...