


Product Datasheet

ACSS2/Acetyl Coa Synthetase Rabbit Monoclonal Antibody (orb548108)

Description	ACSS2/Acetyl Coa Synthetase Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthesized peptide derived from human ACSS2
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	78580 MW
Uniprot ID	Q9NR19
Dilution Range	WB 1:500-1:2000 ICC/IF 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of ACSS2 expression...