


Product Datasheet

PSB9 / LMP2 Rabbit Monoclonal Antibody (orb548121)

Description	PSB9 / LMP2 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, WB
Immunogen	A synthesized peptide derived from human Proteasome 20S LMP2
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	At -20°C for one year. Avoid repeated freezing and thawing
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P28065
Sensitivity	50pg/ml
Dilution Range	WB 1:1000-1:5000ICC/IF 1:50-1:200FC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of Proteasome 20S ...