


Product Datasheet

Calnexin CANX Rabbit Monoclonal Antibody (orb548199)

Description	Calnexin CANX Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human Calnexin
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	67568 MW
Uniprot ID	P27824
Dilution Range	WB 1:1000-1:2000IHC 1:100-1:500ICC/IF 1:50-1:200FC 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of Calnexin expres...