

Product Datasheet

gamma Synuclein Rabbit Monoclonal Antibody (orb548221)

Description gamma Synuclein Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Rat

Conjugation Unconjugated

Tested Applications FC, ICC, IF, IHC, IP, WB

Immunogen A synthesized peptide derived from human gamma Synuclein

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Note For research use only


Isotype Rabbit IgG

Clonality Monoclonal

Uniprot ID **076070**

Dilution Range WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200IP 1:20FC 1:30

Expiration Date 12 months from date of receipt.



Western blot analysis of gamma Synuclein...