


Product Datasheet

NRCAM Monoclonal Antibody (orb548244)

Description	NRCAM Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthesized peptide derived from human NRCAM Cell adhesion, ankyrin-binding protein involved in neuron-neuron adhesion. May play a role in the molecular assembly of the nodes of Ranvier.
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q92823
Sensitivity	10pg/ml
Dilution Range	WB 1:1000-1:5000IP 1:20
Expiration Date	12 months from date of receipt.



Western blot analysis of NRCAM expressio...