


Product Datasheet

AMPK gamma 1 Rabbit Monoclonal Antibody (orb548306)

Description	AMPK gamma 1 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	A synthesized peptide derived from human AMPK gamma 1
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	37579 MW
Uniprot ID	P54619
Dilution Range	WB 1:1000-1:2000IP 1:20FC 1:20
Expiration Date	12 months from date of receipt.



Western blot analysis of AMPK gamma 1 ex...