



Product Datasheet

delta Sarcoglycan SGCD Rabbit Monoclonal Antibody (orb548360)

Description	delta Sarcoglycan SGCD Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	A synthesized peptide derived from human delta Sarcoglycan
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	32071 MW
Uniprot ID	Q92629
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200IP 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of delta Sarcoglyc...



Immunohistochemical analysis of paraffin...