

Product Datasheet

beta I Tubulin TUBB1 Rabbit Monoclonal Antibody (orb548367)

Description beta I Tubulin TUBB1 Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, ICC, IF, IHC, WB

Immunogen A synthesized peptide derived from human beta I Tubulin

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Note For research use only

Isotype Rabbit IgG


Clonality Monoclonal

MW 50327 MW


Uniprot ID [Q9H4B7](#)

Dilution Range WB 1:3000-1:10000 IHC 1:100-1:200 ICC/IF 1:100-1:200 FC 1:100


Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of beta I Tubulin ...



IF analysis of immunocytochemical section...