


Product Datasheet

COX IV COX4I1 Rabbit Monoclonal Antibody (orb548370)

Description	COX IV COX4I1 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, IP, WB
Immunogen	A synthesized peptide derived from human COX IV
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	19577 MW
Uniprot ID	P13073
Dilution Range	WB 1:1000-1:2000IHC 1:50-1:200ICC/IF 1:100-1:500IP 1:20FC 1:20
Expiration Date	12 months from date of receipt.



Western blot analysis of COX IV expressi...