



Product Datasheet

Histone H3.3 H3F3A Rabbit Monoclonal Antibody (orb548426)

Description	Histone H3.3 H3F3A Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human Histone H3.3
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	15328 MW
Uniprot ID	P84243
Dilution Range	WB 1:500-1:1000IHC 1:500-1:1000ICC/IF 1:500-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of Histone H3.3 ex...