

Product Datasheet

Cytokeratin 4 KRT4 Rabbit Monoclonal Antibody (orb548432)

Description Cytokeratin 4 KRT4 Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, IHC, WB

Immunogen A synthesized peptide derived from human Cytokeratin 4

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Note For research use only

Isotype Rabbit IgG

Clonality Monoclonal


MW 57285 MW

Uniprot ID **P19013**

Sensitivity > 5000 cells

Dilution Range WB 1:1000-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of Cytokeratin 4 e...