


Product Datasheet

Histone H3 (mono+di+tri methyl K14) HIST1H3A Rabbit Monoclonal Antibody (orb548492)

Description	Histone H3 (mono+di+tri methyl K14) HIST1H3A Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthesized peptide derived from human Histone H3 (mono+di+tri methyl K14)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	15404 MW
Uniprot ID	P84243, P68431, Q16695
Dilution Range	WB 1:500-1:1000 ICC/IF 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H3 (mon...