

Product Datasheet

Histone H4 (tri methyl K20) HIST1H4A Rabbit Monoclonal Antibody (orb548498)

Description Histone H4 (tri methyl K20) HIST1H4A Rabbit Monoclonal Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen A synthesized peptide derived from human Histone H4 (tri methyl K20)

Form/Appearance Liquid

Concentration Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Storage Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Note For research use only

Isotype Rabbit IgG


Clonality Monoclonal

MW 11367 MW


Uniprot ID **P62805**

Dilution Range WB 1:500-1:2000 IHC 1:50-1:200

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of Histone H4 (tri...