


Product Datasheet

GRK3 Monoclonal Antibody (orb548509)

Description	GRK3 Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human GRK3 Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors.
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P35626
Dilution Range	WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of GRK3 expression...