


Product Datasheet

Histone H2B (yeast) HIST1H2BC Rabbit Monoclonal Antibody (orb548513)

Description	Histone H2B (yeast) HIST1H2BC Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from yeast Histone H2B (yeast)
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	14252 MW
Uniprot ID	P02293
Dilution Range	WB 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H2B exp...