


Product Datasheet

Maltose Binding Protein maleE Rabbit Monoclonal Antibody (orb548516)

Description	Maltose Binding Protein maleE Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from human Maltose Binding Protein
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	43388 MW
Uniprot ID	P0AEX9
Dilution Range	WB 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of MBP expression
...