



Product Datasheet

RFP rfp1 Rabbit Monoclonal Antibody (orb548524)

Description	RFP rfp1 Rabbit Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthesized peptide derived from human RFP
Form/Appearance	Liquid
Concentration	Actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Storage	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	26053 MW
Uniprot ID	Q8ISF8
Dilution Range	WB 1:1000-1:5000 ICC/IF 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of RFP expression ...



Immunofluorescent analysis of 293 cells,...